Optimizing Large-scale Expansion of Mesenchymal Stem Cells in 3D Tide Motion Bioreactors



Geraldine Giap Ying Chiew, Jia Xing Ng, Alan Tin-Lun Lam, Jonathan Han Wei Lim, Shao Liang Lim, Xiangliang Lin

Β.

C.

Source 1

Source 2

Introduction

The future cell therapy demands high quantities of Mesenchymal stem cells (MSCs) ranging from 10 million to more than 200 million cells per dosage. Conventional expansion of MSCs on plastic wares become impractical when large dosages of more than 50 million cells are required. The use of bioreactors which combines scaling-up ability, process control and automation is the primary solution for this need. Many bioreactors are facing issues in supporting MSC cultures due to complications in balancing the need for proper mixing of media with the need to minimize shear stress as well as the inability to separate cells from micro/macro carriers with high cell yield and viability.

Esco Aster has leveraged on the Tide Motion bioreactors to establish a robust and scalable platform using macrocarriers to meet the demands for future clinical therapies. MSCs isolated from different tissues sources were seeded and allowed to expand within PET macrocarriers. Throughout culture periods, cell culture conditions were monitored, with bioprocess parameters such as glucose consumption and pH levels measured to ensure proper scale-up. Key issues such as cell seeding densities, media culturing conditions and improved bioprocess parameters needed for optimal stem cell systems were studied in our system. Overall, we present our process optimization with quality controls and release criteria of functional and phenotypic characteristics for the translation of academic/industrial R&D into bench sale for future clinical trials and commercialization process.



Tissue Origin	Cells Seeded	Cells Harvested	Fold Increase	Culture Days
UC-MSC	1.96E+07	2.41E+08	12.33	7
WJ-MSC	2.55E+07	2.59E+08	10.14	6
CL-MSC	2.85E+07	1.89E+08	6.64	6
BM-MSC	1.70E+07	1.66E+08	9.79	7
AD-MSC	2.50E+07	2.30E+08	9.15	9

(A) Lower attachment efficiency on the carriers was observed when UC-MSCs were seeded in hPL supplemented media. (B) Despite the lower attachment efficiency observed, hPL supplemented media gave rise to the greatest expansion of cells across MSCs from different tissue origins.

(A) Successful expansion of MSCs from different tissue origins have been demonstrated in the system. Common tissue origins of MSCs such as umbilical cord-derived MSC (UC-MSC), Wharton Jelly-derived MSC (WJ-MSC), bone marrow-derived MSC (BM-MSC) and adiposederived MSC (AD-MSC) generally show a 10-fold expansion within 7 days of culture. The table above lists the recommended seeding density for the various tissue sources of MSC optimized using commercial sources of the cells. (B) Different sources of MSCs from the same tissue origin were tested for their ability to expand from a low seeding density. Florescence imaging of the cells on the carriers at the end of the culture period show successful expansion to a confluent level by BM-MSCs obtained from source 2 unlike MSCs from source 1. As differences in growth parameters were observed between MSCs from different sources, optimization of parameters such as seeding density and media change regime is usually conducted for each source of MSC prior to expansion in bioreactors of a larger scale.

(C) UC-MSCs were adapted to the different media types over 2 passages before seeding onto BioNOC II carriers for monitoring of proliferation rates. (D) Fluorescent images show the expansion of the cells on the carriers over the culture period. Green: Fluorescein diacetate (cytoplasm of live cells), Blue: Hoechst 33342 (nucleus), Red: propidium iodide (dead cells)

Optimization of Harvesting Process

Testing of Different Enzymes

	Accumax	Collagenase	TryPLE Express	Trypsin
Harvesting (%)	87	68.3	78.3	56.3
Viability (%)	95.6	73.2	95.8	88.9

Mass Harvesting from CelCradle[™]

mass harvesting of UC-MSCs from BioNOC II carriers. Fluorescein diacetate (FDA) was used to stain for remaining live cells on BioNOC II carriers were sampled throughout the harvesting process.

Harvest	Live Cells Harvested	Viability	Total Live Cells Harvested
Enzyme 1	14,625,000	94.35%	
Enzyme 2 + Wash 1	181,500,000	96.80%	233,375,000
Wash 2 – 5	37,250,000	95.82%	

Quality Checks

Conc	lusion

Classification	Parameters	Procedure	Reference value
Quantity	Cell Count	Automated cell counter	Depending on dose per vial (Range from 10 million to 1 billion)
Viability	Cell viability	Automated cell counter	>90%
Identity	Karotype	Cytogeneticists study	Normal Karyotype
Identity	Morphology	Visual microscopy observations	Fibroblastic spindle shape
Identity	Trilineage differentiation	Differentiation assay (Adipogenic, Chrondrogenic, Osteogenic)	Pos
Identity/ Purity	Stemness	Immune phenotyping by FACS	≥ 95% CD105, CD73, CD90 ≤ 2% CD34, 45, HLA DR, CD11b, CD19
Safety	Sterility <usp71></usp71>	Automated microbial detection system	Neg
Safety	Mycoplasma <usp63></usp63>	PCR	Neg
Safety	Endotoxin <usp85></usp85>	Limulus Amebocyte Lysate	<0.25EU/kg
Safety	Viral (HIV, HBV, HCV, HPV)	Viral protein antigen (ELISA)	Neg
Safety	Tumorigenicity tests	Telomerase enzyme activity	RTA<1.2
Potency	Cytokine protein expression assay	eg IL-6 protein secretion (ELISA)	
Potency	Cytokine gene expression assay	eg TNF- α gene expression (q-PCR)	
Potency	MSC migration/adhesion	Scratch test/ Transwell assay	
Potency	Immunomodulatory assay	IDO expression and activity assay	

 Tide Motion bioreactors 	Factor/Attribute	Performance of MSCs	
consists of a series of single- use, loose packed-bed	Cell Types	Suitable for various tissue origins of MSCs (UC, WJ, BM and AD)*	
 They serve as a platform for high capacity manufacturing of 	Cell Attachment	High seeding efficiency >90% Fibronectin coating for serum free media (suitable for cGMP production)	
MSCs in single batches.	Cell Growth & Monitoring	Ease of visualizing cells on matrix via dye stain Confluence achieved at day 6-8 with ~10X expansion	
Small scale systems support early developmental stages of	Cell Harvest and Media Collection	>90% cells harvested with suitable enzymes Permits collection of media for processing of exosomes or conditioned media	
the bioconversion process, allowing for optimization of	Quality of cells	Stemness and trilineage differentiation of MSCs preserved Viability of harvest >90%, healthy cells obtained after harvest	
parameters at low costs.			

Current seeding at 2e7 cells per CelCradle[™] to obtain 2e8 cells

Projected to obtain 4e9 cells in 2L TideXcell[™] and 6e11 cells in 300L TideXcell^{™**} *Growth parameters of other tissue origins of MSCs scheduled for development **Final density will vary based on age, source of stem cells and media type used

© 2019 Esco Aster Pte Ltd. All rights reserved. All trademarks are the property of Esco unless otherwise stated.

ESCO ASTER CDMO SERVICES & TIDE MOTION BIOREACTORS

Find out more at http://www.vaccixcell.com or http://www.escoaster.com For enquiries, contact Esco Vaccixcell (Bioprocessing Tools): mail@vaccixcell.com or Esco Aster (cGMP CDMO): mail@escoaster.com