

Scalable Production of Modified HEK293T Extracellular Vesicles Using Adherent Packed Bed Bioreactor Cell Culture System




Choon Keong Lee, Claudine Ming Hui Lim, Winnie Faustinelie, Desy Silviana, Xiangliang Lin


Summary

Our study highlights the efficacy of our proprietary adherent packed bed bioreactor for a scalable production of extracellular vesicles (EVs). Leveraging on the innovative ‘Tide Motion’ principle, which alternately exposes the cells to air and nutrients to maintain the cell culture system. In a commercial endeavor, we engineered HEK293T cells to express a specific protein on the surface. The cells were seeded in our smallest bioreactor, the CelXrocker®, as a proof-of-concept evaluation. We successfully demonstrated the feasibility of cultivating the cells in our bioreactor and the CelXrocker® may minimize the cost and time for early-stage process development. Subsequently, the cell culture were scale up in the CelCradle® and to support the high-density cell growth, a perfusion system was attached to the bioreactor (CelCradle X®). Our data shows that the bioreactor can support the high-density cell growth of the modified HEK293T cells, and this significantly increased the EVs production yield. A correlation between the cell culture parameters and the quality of EVs can be established. We carried out downstream processing to purify the EVs and analytical characterizations were carried out to confirm the identity of the EVs.

General Process Development Steps



Upstream Process Development



Downstream Process Development and Analytical Characterization

